

Comparing Timing Resolution Improvements using the Tapered-Sinc Interpolation Algorithm for Various Readout Methods



2024 MTV Workshop

March 27

Kim Meagher

Undergraduate Student, University of Michigan

T. E. Maurer¹, S. D. Clarke¹, S. A. Pozzi¹

¹University of Michigan



Introduction and Motivation

- Detection imaging systems used for nonproliferation purposes require high precision timing
- *NNSA mission relevance and other applications*
 - Improving timing resolution will improve...
 - **Imaging systems that can locate orphan sources for nonproliferation and safeguards purposes**
 - Positron Emission Tomography (PET) systems to evaluate biochemical processes
- Digitizers for data acquisition have fixed sampling rates
 - Can limit coincidence timing resolution (CTR) of fast scintillators like organic scintillators




Orphan Cs-137 source located via radiation imaging systems in Australian Outback

Photo credit: Radford, Antoinette. "How a Tiny Radioactive Capsule Was Found in Australia's Vast Outback." *BBC News*, BBC, 1 Feb. 2023.

T-Sinc Function Motivation

- Post processing interpolation algorithms can increase sampling rate of pulse, ultimately improving coincidence timing resolution (CTR)
- Fill in rising edge of pulse
- Use constant fraction discriminator (CFD) to measure CTR
 - Improve CTR via interpolation and increased sampling rate with t-sinc algorithm
- Investigate functionality on Photomultiplier Tube (PMT) vs. Silicon Photomultiplier (SiPM) readout methods



Fast pulse interpolated using t-sinc algorithm

Photo Credit: W. K. Warburton and W. Hennig, "New Algorithms for Improved Digital Pulse Arrival Timing With Sub-GSps ADCs," in IEEE Transactions on Nuclear Science, vol. 64, no. 12, pp. 2938-2950, Dec. 2017, doi: 10.1109/TNS.2017.2766074.

T-Sinc Function Background

T-Sinc is the tapered sinc function:

- $\text{sinc}(x) = \sin(x)/x$
- $\text{tsinc}(i) = \text{sinc}(i\pi/N)\exp(-(i/T)^2)$
- $$g(j, k) = \sum_{i=0}^{L-1} (y(j-i)\text{tsinc}(iN+k) + y(j+1+i)\text{tsinc}((i+1)N-k))$$

N - Number of interpolated points

T - Tapering constant

L - Number of lobes

W. K. Warburton and W. Hennig, "New Algorithms for Improved Digital Pulse Arrival Timing With Sub-GSps ADCs," in IEEE Transactions on Nuclear Science, vol. 64, no. 12, pp. 2938-2950, Dec. 2017, doi: 10.1109/TNS.2017.2766074.



Coincidence Timing Resolution (CTR) Background

- CTR quantifies system resolution
 - **Readout method**
 - Photomultiplier Tube
 - Silicon Photomultiplier
- Annihilation photons from centered ^{22}Na source should theoretically hit detectors at the same time/consistent difference in time
 - Deviance implies non-zero timing resolution
 - Quantified using full-width at half maximum




Photo Credit: Pots, Rosalinde & Gundacker, Stefan & Auffray, Etienne. (2020). "Exploiting Cross-Luminescence in BaF₂ for Ultrafast Timing Applications Using Deep-Ultraviolet Sensitive HPK Silicon Photomultipliers." *Frontiers in Physics*. 8. 10.3389/fphy.2020.592875.

Coincidence Timing Resolution Experiment Set-Up

Two experiments performed comparing different readout methods...

1. Photomultiplier tubes:

- Bare 6 mm edge organic glass scintillator (OGS)* cubes tested and coupled to 38.1 mm diameter Hamamatsu H3178-51 PMTs

2. Silicon photomultipliers:


- Teflon wrapped 6 mm edge OGS* cubes tested
- SiPM readout on single board

Both experiments used...

- Cs-137 source for calibration
- 2.5 GHz oscilloscope for acquisition
- Na-22 centered between detectors to measure annihilation photons



PMT setup



SiPM setup

*  Sandia National Laboratories



T-Sinc Analysis Procedure

1. Using the CoMPASS data acquisition system to acquire fast pulses from 2.5 GHz oscilloscope, find the CFD time difference between Channel 1 and Channel 2 for every set of coincidence pulses.
 - a. CFD fraction = 20%
 - b. Threshold set at 50mV
2. Histogram the time difference
3. Repeat the above steps with down-sampled data and interpolated data.




T-Sinc Analysis Procedure

Original Sampled Pulse
Sampling Frequency = 2.5 GHz
Downsampling factor $N=3$




Example Reconstructed Pulses

Readout: PMT




Readout: SiPM




T-Sinc Function on PMT Data $N=5$

Actual Pulse




FWHM = 0.660 ns

Down-Sampled Pulse



FWHM = 0.848 ns


Interpolated Pulse



FWHM = 0.688 ns


T-Sinc Function on SiPM Data $N=3$

Actual Pulse




FWHM = 0.362 ns

Down-Sampled Pulse



FWHM = 0.324 ns

Interpolated Pulse



FWHM = 0.374 ns

Results Summarized

	3.81 cm Diameter PMT	Single Board SiPM
Actual pulse FWHM (ground truth) (ns)	0.66	0.36
Down-sampled pulse FWHM (ns)	0.85	0.32
Interpolated pulse FWHM (ns)	0.69	0.37
Actual and interpolated FWHM percent difference (%)	4.5	2.8
Down-sampled and interpolated FWHM percent difference (%)	23	1.6



Conclusion and Impacts

- T-sinc improves coincidence timing resolution by 20% on photomultiplier tube readout with organic glass cubes
 - Does not have the same impact on silicon photomultiplier data
 - Perhaps due to differing pulse shape of SiPM pulses
- The t-sinc algorithm gives the potential of a **higher sampling frequency** while using low frequency readout methods
 - Impacts experiments that require high precision timing such as coincidence timing resolution
- Improved sampling rate in post processing can improve **radiation imaging systems** for **safeguarding purposes** in line with NNSA mission



Next Steps

- Investigate t-sinc parameters
 - N - Number of interpolated points
 - T - Tapering constant
 - L - Number of lobes
- Investigate other post-processing methods
 - Step function
- Investigate other equipment
 - Various diameter PMTs
 - Various data acquisition systems



Acknowledgements



The Consortium for Monitoring, Technology, and Verification would like to thank the DOE-NNSA for the continued support of these research activities.



This work was funded by the Consortium for Monitoring, Technology, and Verification under Department of Energy National Nuclear Security Administration award number DE-NA0003920

